



fittleworth

# Laryngectomy and Tracheostomy

Catalogue



## Opening Hours:

8am to 8pm Mon to Fri

9am to 1pm Sat

**National:** 0800 378 846

**Scotland:** 0800 783 7148

**Nurse Line:** 0800 378 413

**Text:** 07800 005 658

**Online:** [www.fittleworth.com](http://www.fittleworth.com)

Issue 12

# Contents

## SECTION 1: Heat and Moisture Exchangers (HMEs) & Baseplates

Blom-Singer®	4 - 7
FAHL	8 - 11
Provox®	12 - 13
SpiroTect™	15 - 18
Phoniq	19

## SECTION 2: Tracheostomy and Laryngectomy Accessories

Insight Medical Products	21 - 23
Adhesives and Adhesive Removers	24
Barrier Creams and Lubricants	25
Barrier Film Wipes	26
Voice Prosthesis Brushes	27
Tube Cleaning Brushes	28
Dressings	29 - 31
Protectors	32 - 34
Shower and Foam Protectors	35
Trachi-Naze +Plus Nasal Restoration System	36
Tube and Collar Holders	37 - 39
Tube holders	40
Other Accessories	41 - 43

## SECTION 3: For Healthcare Professionals


FAHL Laryngectomy Tubes	45
Kapitex Trachi-Naze® Plus Tubes	46
SpiroTect™ Laryngectomy Tubes	47
Shiley™ Tracheostomy Tubes	48 - 50
Provox Larytube	51
Therabite® Jaw Motion Rehab System™	51

## SECTION 4: About Home Delivery

How does this service work?	53
Reordering list	54 - 55

For a complete and up-to-date list of products available, please refer to the Drug Tariff.

All text, images and design in this publication are subject to copyright.



# Heat and Moisture Exchangers (HMEs) & Baseplates



## Blom-Singer® HME System

### Blom-Singer® ClassicFlow™ HME System Starter Kit

Convenient seven-day starter kit for those laryngectomees wanting to try the ClassicFlow HME System (now also includes the TruSeal Contour Low Profile Adhesive Housings).

Order ref: BE1055EZ

Pack size: 1



BE1055EZ

### Blom-Singer® HumidiFilter® Holder

This lightweight, reusable housing attaches to the front of the tracheostoma adhesive housing. Use in conjunction with HumidiFilter® Holder Foam Inserts (BE 1080) to filter, heat and moisten inspired air.

Order ref: BE1060

Pack size: 1



BE1060

### Blom-Singer® Replacement Foam Filters (for non-ATSV II users)

Specially treated foam discs that fit ONLY with the HumidiFilter Holder (BE1060). These filters DO NOT fit with the ATSV II HumidiFilter Cap (BE 1010).

Order ref: BE1080

Pack size: 30



BE1080

## Blom-Singer® Hands-Free Speech

### Blom-Singer® Hands-Free (ATSV II) System

The assembled two-piece valve inserts into the tracheostoma adhesive housing that is applied to the skin surrounding the stoma. Consists of a removeable diaphragm/faceplate and body.

Order ref: BE8025H

Pack size: 1



### Blom-Singer® HumidiFilter® ATSV II Cap and 7 Foam Filters

The HumidiFilter® ATSV II Cap is a twist-on plastic filter cap that fits over the ATSV II. It is designed to help keep the ATSV II Foam Filters in place and come with 7 specially treated HME foam Filters.

Order ref: BE1010

Pack size: 1



### Blom-Singer® Replacement Foam Filters

Specially treated foam discs that fit ONLY with the HumidiFilter® ATSV II Cap (BE1010).

Order ref: BE1030

Pack size: 30



### Blom-Singer® Replacement Diaphragm/Faceplate for ATSV II

Replacement diaphragm/faceplate for Blom-Singer® ATSV II Valve.

Order ref: BE8026H

Pack size: 1



## Blom-Singer® Adhesive Baseplates

### Blom-Singer® TruSeal® Contour™ Adhesive Housings

The TruSeal Contour Adhesive Housing is a lightweight, flexible material offered in three shapes that fit various peristomal anatomies. The TruSeal Contour is a disposable, all-inclusive adhesive tape-disc/housing that eliminates the need for a tracheostoma valve housing, tape or adhesive foam disc. It is designed to hold the ATSV II, HME Cartridge or shower guard over the stoma.

#### Round

Order ref: Regular	BE6070	<input type="checkbox"/>
Pack size: 30		

Order ref: Low Profile	BE6074	<input type="checkbox"/>
Pack size: 30		



#### Oval

Order ref: Regular	BE6071	<input type="checkbox"/>
Pack size: 30		


Order ref: Low Profile	BE6075	<input type="checkbox"/>
Pack size: 30		



#### Standard

Order ref: Regular	BE6072	<input type="checkbox"/>
Pack size: 30		

Order ref: Low Profile	BE6076	<input type="checkbox"/>
Pack size: 30		



### Blom-Singer® TruSeal® Adhesive Housing

The TruSeal Adhesive Housing is a lightweight, disposable, all inclusive adhesive tape-disc/housing that eliminates the need for a tracheostoma valve housing, tape or adhesive foam disc. It is designed to hold the ATSV II HME cartridge or shower guard over the stoma.

Order ref: BE6054	<input type="checkbox"/>
Pack size: 30	

## Blom-Singer® Tracheostoma Housing and Adhesive Discs

### Blom-Singer® Tracheostoma Valve Housing

The Tracheostoma Valve Housing is a reusable, flexible housing that attaches the ATSV II, HME cartridge or shower guard over the stoma. It is applied with separate adhesive tape discs and silicone adhesive. It is available in two sizes, standard and large.



Order ref: Standard	BE6038	<input type="checkbox"/>
Pack size: 1		

Order ref: Large	BE6039	<input type="checkbox"/>
Pack size: 1		

### Blom-Singer® Adhesive Discs

These double-sided foam discs provide the optimal attachment of the tracheostoma valve housing to your stoma. They are available in two different sizes and are supplied in packs of 30 discs.

#### Heavy Duty Adhesive Tape Discs

Order ref: Standard	BE6034	<input type="checkbox"/>
Pack size: 30		

Order ref: Large	BE6035	<input type="checkbox"/>
Pack size: 30		

#### Thin Adhesive Tape Discs

Order ref: Standard	BE6043	<input type="checkbox"/>
Pack size: 30		

Order ref: Large	BE6044	<input type="checkbox"/>
Pack size: 30		

#### Regular Adhesive Tape Discs

Order ref: Standard	BE6041	<input type="checkbox"/>
Pack size: 30		

Order ref: Large	BE6042	<input type="checkbox"/>
Pack size: 30		

## FAHL HME Systems

### Laryvox® Pro HME and Laryvox® Pro HighFlow HME

#### Features and Benefits

##### Tinted Top

- The Laryvox® Pro HME and Laryvox® Pro HME HighFlow have tinted tops to reduce visibility
- The tinted top has a low glare surface finish to reduce light reflections

##### 24 Hour Use

- The Laryvox® Pro HME and Laryvox® Pro HME HighFlow are daily disposable, 24 hour use devices
- Side vents aid night-time use

##### Speech Button

- One touch seals the valve to speak
- Release the valve to return to breathing

#### Laryvox® Pro HME

- The Laryvox® Pro HME provides normal airflow resistance
- Applicable for all patients following laryngectomy surgery
- Used for normal physical exertion

Order ref: 49841

Pack size: 30



#### Laryvox® Pro HME HighFlow

- Low airflow resistance therefore easy to breathe through
- Applicable for all patients following laryngectomy
- Ideal starting point for all new users of HMEs or those who are used to wearing bibs

Order ref: 49840

Pack size: 30



FAHL

Andreas Fahl®  
Laryvox® Pro HME



FAHL

Andreas Fahl®  
Laryvox® Pro HME HighFlow

*With friendly permission of Andreas Fahl Medizintechnik-Vertrieb GmbH, Colgogne, Germany.*

National: 0800 378 846 | Scotland: 0800 783 7148 | Nurse Line: 0800 378 413 | Text: 07800 005 658



## FAHL HME Systems

### Laryvox® HME

An excellent HME cassette, designed to reduce mucus production, coughing and chest infections. 24 hour use also helps to improve sleep and reduce levels of fatigue in patients post-laryngectomy.



Order ref: 49800  
Pack size: 30



### Laryvox® O2 HME

Laryvox® O2 HME provides a unique combination of Laryvox® HME with an O2 connection, ideal for the patient who requires supplementary oxygen, for example in the post-operative setting. Designed to reduce mucus production.



Order ref: 49802  
Pack size: 30



### Laryvox® HME HighFlow

The Laryvox® HME HighFlow offers low airflow resistance so is very easy to breathe through. This makes it the perfect choice for individuals who are new users of HME, in particular those who may be switching from a foam stoma protector or a bib.



Order ref: 49810  
Pack size: 30



*With friendly permission of Andreas Fahl Medizintechnik-Vertrieb GmbH, Colgogne, Germany.*

National: **0800 378 846** | Scotland: **0800 783 7148** | Nurse Line: **0800 378 413** | Text: **07800 005 658**

## FAHL Adhesive Baseplates

### Laryvox® Tape Hydrosoft

The hydrosoft baseplate is made of a soft, hydrocolloid material and incorporates a very mild adhesive.

- Particularly skin friendly
- Comfortable on application
- Ability to enclose fluids/secretions, ideal for post-operative application
- Not suitable for use in combination with shower guard

Order ref: Round	48130	<input type="checkbox"/>
Pack size: 15		
Order ref: Oval	48230	<input type="checkbox"/>
Pack size: 15		



48130



48230

### Laryvox® Tape Hypoallergen

A unique comfortable and flexible baseplate. Designed for patients with skin sensitivities or skin reactions but also offers strong adhesion.

- Adapts to movement
- Available in three sizes

Order ref: Round	48050	<input type="checkbox"/>
Pack size: 15		
Order ref: Oval	48150	<input type="checkbox"/>
Pack size: 15		
Order ref: XL Oval	48250	<input type="checkbox"/>
Pack size: 15		



48050

*With friendly permission of Andreas Fahl Medizintechnik-Vertrieb GmbH, Colgogne, Germany.*

National: **0800 378 846** | Scotland: **0800 783 7148** | Nurse Line: **0800 378 413** | Text: **07800 005 658**

## FAHL Adhesive Baseplates

### Laryvox® Tape Comfort

A comfortable, tinted base plate which adapts well to the contours of the skin. With the strongest level of adhesion it makes an ideal choice for patients experiencing issues with baseplate retention or who wish to trial a Hands-Free system e.g. the Blom-Singer ATSV II.



Order ref: Oval	48240	<input type="checkbox"/>
Pack size: 15		

Order ref: XL Oval	48340	<input type="checkbox"/>
Pack size: 15		

### Laryvox® Tape Extra Fine

A unique transparent, strong, waterproof baseplate that provides long lasting adhesion. The ultra thin, transparent material adapts well to all tracheostomal conditions, making it an ideal choice for those patients with an irregular or recessed stoma.



Order ref: 48400	<input type="checkbox"/>
Pack size: 15	

### Laryvox® Tape Standard and Flexible (Round, Oval and XL Oval)

- Transparent appearance with water-repellent effect
- Non-invasive
- High adhesive strength on the skin around the tracheostoma
- Suitable for problematic tracheostoma situations
- Suitable for normal physical activity
- No friction inside the tracheostoma
- With 22mm housing to be combined with any HME and speaking valve with 22 mm adapter (KOMBI)
- Can be combined with special silicone tubes (LARYNGOTEC® KOMBI CLIP or LARYNGOTEC® KOMBI CLIP LINGO) if additional tracheal stabilisation is needed

#### - Standard

Order ref: Round	48100	<input type="checkbox"/>
Pack size: 15		

Order ref: Oval	48200	<input type="checkbox"/>
Pack size: 15		

Order ref: XL Oval	48300	<input type="checkbox"/>
Pack size: 15		

#### - Flexible

Order ref: Round	48120	<input type="checkbox"/>
Pack size: 15		

Order ref: Oval	48220	<input type="checkbox"/>
Pack size: 15		

With friendly permission of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne, Germany.

National: 0800 378 846 | Scotland: 0800 783 7148 | Nurse Line: 0800 378 413 | Text: 07800 005 658

## Provox® Freehands FlexiVoice HME & Provox® Luna HME System

The FreeHands FlexiVoice offers both freedom to speak hands-free and with finger occlusion when wanted. For hands-free speech you can choose between 3 different strengths; light, medium and strong. The choice depends on how much airflow and pressure you need for speaking and breathing.

- Speak hands-free speech or with finger occlusion
- Use different devices for different situations
- Lock for physical activity
- Fast set-up
- Reduces mucus production and coughing

### Provox® FreeHands FlexiVoice - Light

Order ref: 8161   
Pack size: 1

### Provox® FreeHands FlexiVoice - Medium

Order ref: 8162   
Pack size: 1

### Provox® FreeHands FlexiVoice - Strong

Order ref: 8163   
Pack size: 1

### Provox® FreeHands FlexiVoice HMEs - Moist

Order ref: 8220   
Pack size: 1

### Provox® FreeHands FlexiVoice HMEs - Flow

Order ref: 8221   
Pack size: 1

### Provox® Luna HME

Provox Luna is a new night-time solution that helps you sleep comfortably, soothe your skin and improve your lung health.

- Soft and smooth for night-time comfort
- Made of hydrogel to soothe and cool the skin
- Superior humidification compared to other HMEs
- Low breathing resistance for easy breathing at night



8013

Order ref: 8013   
Pack size: 30

## Provox® HME System

### HME Filter Cassette

Provox® HMEs reduce coughing and mucus production resulting in best possible humidification and comfortable breathing. Can be used day or night.

### XtraMoist HME

Helps to keep your lungs healthy by humidifying and filtering the air that you breathe. Wear it at night so you can get the best possible humidification and breathe comfortably.

Order ref: 7290  
Pack size: 30



### XtraFlow HME

Helps to keep your lungs healthy by humidifying and filtering the air that you breathe. Put it on during the day when you are physically active.

Order ref: 7291  
Pack size: 30



### Micron® HME

The Micron® HME gives you the best filtering possible. Use it in crowded places like schools, supermarkets or on public transport, especially during the flu season.

It can also be helpful whilst gardening, particularly during the high pollen season.

Order ref: 7248  
Pack size: 30



### **StabiliBase™**

Provox® StabiliBase™ is a strong adhesive that gives you better stability, especially if your stoma is deep.

You can use it during the day when you speak a lot, or speak hands-free with Provox® FreeHands FlexiVoice.

### **Flexiderm Baseplate**

Designed to hold the Provox® HME Cassette, Provox® FreeHands HME and Provox® ShowerAid in front of the tracheostoma.

The baseplate provides an optimised seal to the skin around the stoma.

There are 4 different types of adhesives, with 3 different adhesive properties and two different shapes.

### **Provox® XtraBase**

- For speaking with finger occlusion or hands-free
- For a flat or deeper stoma
- Adheres well to the skin

### **Provox® LaryClip**

The Provox® LaryClip is an alternative to the Provox® TubeHolder. It is a two-piece system that helps to optimize the air-tight attachment of Provox® LaryButton and Provox® LaryTube to the stoma.

### **Provox® OptiDerm**

Provox® OptiDerm is a soft adhesive that gives you a gentler seal around your stoma. Use it during the night to give your skin some relief or during the day if your skin feels sore or irritated.

### **Provox® Regular**

- For oesophageal speech
- For speaking with electrolarynx
- Easy to remove from the skin

## Provox® Adhesive Baseplates

### Plus Size Baseplates

*Stronger adhesion for difficult stomas*

Provox® Adhesive Baseplates are designed to hold the Provox® HME Cassette, Provox® Freehands HME and Provox® Shower Aid in front of the tracheostoma, and to provide an optimized seal to the skin around the stoma.

There are four different types of adhesives, with three different adhesive properties and two different shapes (oval and round).

#### Provox® XtraBase - Convex

Order ref: 7265   
Pack size: 20

#### Provox® StabiliBase - Convex

Order ref: 7289   
Pack size: 15

#### Provox® OptiDerm Adhesive - Plus

Order ref: 7332   
Pack size: 20

#### Provox® OptiDerm Adhesive - Round

Order ref: 7255   
Pack size: 20

#### Provox® OptiDerm Adhesive - Oval

Order ref: 7256   
Pack size: 20

#### Provox® FlexiDerm Adhesive - Plus

Order ref: 7331   
Pack size: 20

#### Provox® FlexiDerm Adhesive - Round

Order ref: 7253   
Pack size: 20

#### Provox® FlexiDerm Adhesive - Oval

Order ref: 7254   
Pack size: 20

## SpiroTect™ Systems

All SpiroTect™ HME cassettes contain antibacterial coated foam which can help reduce infection and filter inhaled air. All HME cassettes have a fully tinted design, are small, discreet, easy to use and interchangeable with most other devices on the market.

### SpiroTect™ Touch2Speak Cassette, Antibacterial

- Suitable for day and night use
- Occlusion is easily obtained by light finger touch
- Ideal for post-operative Laryngectomees



PTC1004

Order ref: PTC1004  
Pack size: 30



### SpiroTect™ Push2Speak Cassette, Antibacterial

- Suitable for day and night use
- Occlusion is easily obtained by pressing the top of the cassette
- Ideal for less dexterous Laryngectomees



PTC1005

Order ref: PTC1005  
Pack size: 30



### SpiroTect™ Wear2Protect Cassette, Antibacterial

- Ideal for night use
- For use by Laryngectomees without speech ability
- Perfect for Laryngectomees with oesophageal speech
- Occlusion cannot be obtained



PTC1006

Order ref: PTC1006  
Pack size: 30





## SpiroTect™ Systems

All SpiroTect™ Baseplates are hypoallergenic, lightweight, provide excellent adhesion, prevent air leakage and make speech easier.

### SpiroTect™ Secure Convex Baseplate

- Unique honeycomb convex base for stronger yet flexible overall structure
- Developed for more precise and secure positioning
- Prevents air leakage with enhanced and supportive seal
- Showerproof



PTW700

Order ref: PTW700



Pack size: 20

### SpiroTect™ Fabric Baseplate with Tube Connectors

- Round - PTW103

The SpiroTect™ Fabric Baseplate with Tube Connectors Round (PTW103) are skin friendly fabric baseplates that secure Laryngectomy Tubes at the stoma site.

They are specifically designed to provide comfort, security and prevent air leakage.

- Soft, flexible material
- Flexible centre for secure positioning of laryngectomy tube
- Flexible connectors to attach the laryngectomy tubes
- Can be used for up to 50 hours



PTW103

Order ref: PTW103



Pack size: 20

## SpiroTect™ Systems

### SpiroTect™ Fabric Baseplate

- Oval - PTW100
- Round - PTW101

- Soft fabric material with soft cuff
- Can be used for up to 50 hours
- Showerproof



Order ref: PTW100	<input type="checkbox"/>
Pack size: 20	
Order ref: PTW101	<input type="checkbox"/>
Pack size: 20	

### SpiroTect™ Clear and Strong Baseplate

- Oval - PTW300
- Round - PTW301

- Discreet, transparent design
- Lightweight and flexible
- Showerproof



Order ref: PTW300	<input type="checkbox"/>
Pack size: 20	
Order ref: PTW301	<input type="checkbox"/>
Pack size: 20	

### SpiroTect™ Sensitive Baseplate

- Oval - PTW200
- Round - PTW201


- Soft hydrocolloid material
- Adhesive free
- Ideal for sensitive skin and post-op use



Order ref: PTW200	<input type="checkbox"/>
Pack size: 20	
Order ref: PTW201	<input type="checkbox"/>
Pack size: 20	

### SpiroTect™ Extra Sticky Baseplate

- Soft, foam material
- Very flexible
- Showerproof



Order ref: PTW501	<input type="checkbox"/>
Pack size: 20	

## Phoniq Systems

Code:	Product	Pack Size:
LAPBP0001	Phoniq Baseplate HME System Secure Oval	20 <input type="checkbox"/>
LAPBP0002	Phoniq Baseplate HME System Secure Round	20 <input type="checkbox"/>
LAPBP0003	Phoniq Baseplate HME System Secure Hydrocolloid Oval	20 <input type="checkbox"/>
LAPBP0004	Phoniq Baseplate HME System Secure Hydrocolloid Round	20 <input type="checkbox"/>
LAPFT0001	Phoniq DigiClose Comfort HME Filter	30 <input type="checkbox"/>
LAPFT0002	Phoniq DigiClose Active HME Filter	30 <input type="checkbox"/>
LAPFT0003	Phoniq DigiClose Comfort HME Filter Beige	30 <input type="checkbox"/>
LAPFT0004	Phoniq DigiClose Active HME Filter Beige	30 <input type="checkbox"/>



# Tracheostomy and Laryngectomy Accessories



## Insight Medical Products

### Tracheostomy Tube Holder

- Large - TH100
- Child - TH150



The Insight Tracheostomy Tube Holder uses the highest quality materials to provide a safe and secure adhesion for the tracheostomy tube.


- Adult and Paediatric sizes
- Quality hook/loop fastening
- Cough reflex section, allowing movement and comfort

Order ref: TH100   
Pack size: 10

Order ref: TH150   
Pack size: 10

### Advadraw (T) Dressing

- 80 x 80mm



- Cushioning effect on wound
- Low adherence
- Reduced trauma on removal
- Conformable comfortable protection
- Suitable for low/medium Levels of exudate

Order ref: CR4416   
Pack size: 20

### Advazorb Dressings

- Advazorb Plus  
CR/3900 - 80 x 80mm
- Advazorb Light  
CR/4212 - 80 x 80mm



Made from a lightweight, soft and conforming low-adherent foam with high fluid absorbance and retention characteristics.

Order ref: CR/3900   
Pack size: 20

Order ref: CR/4212   
Pack size: 20

### Tracheostomy Swabs

- 6mm - TCS/950
- 9mm - TCS/900



The Insight Cleaning Swabs are perfect for cleaning all types of PVC and silicone tracheostomy or laryngectomy tubes.

Order ref: TCS/950   
Pack size: 10

Order ref: TCS/900   
Pack size: 10

## Insight Medical Products

### Sofshield Laryngectomy and Tracheostomy Drawstring Protectors

Sofshield Laryngectomy Protectors can be washed a minimum of 10 times, whilst still retaining filtration efficiency, offering a significant cost saving over alternative protectors that still use foam inserts, also eliminating the need for costly additional cosmetic covers.

All Sofshield Protectors have an antibacterial coating to help prevent infections.



Order ref: AS3698 - Small - White   
Pack size: 10

Order ref: AS3777 - Small - Blue   
Pack size: 10

Order ref: AS3778 - Small - Beige   
Pack size: 10

Order ref: AS4056 - Medium - Blue   
Pack size: 10

Order ref: AS4051 - Medium - Black   
Pack size: 10

Order ref: AS3669 - Large - White   
Pack size: 10

Order ref: AS3764 - Large - Blue   
Pack size: 10

Order ref: AS3765 - Large - Beige   
Pack size: 10

### Sofnex Laryngectomy Shield

The Sofnex laryngectomy shield is ideal for wearing with a cravat. All Sofnex Protectors have an antibacterial coating to help prevent infection. Sofnex can be washed at least 10 times with no loss of filtration efficiency or distortion to the product.



Order ref: CR3856 - White   
Pack size: 10

Order ref: CR3859 - Blue   
Pack size: 10

Order ref: CR3860 - Beige   
Pack size: 10

### Insil Adhesive Tape

This versatile tape is backed with a silicone adhesive. When tape is removed, it is extremely gentle on the skin.



Order ref: INSIL/200 - 2cm x 3m   
Pack size: 1

## Insight Medical Products

### Cascade Shower Protector


Cascade is a soft waterproof PVC protector that provides a barrier to prevent water from entering the stoma, thereby improving the patient's confidence to shower routinely.



Order ref: AS/3835   
Pack size: 10 (5 packs of 2)

### Purifoam Laryngectomy Filters

Purifoam Filters use a high quality medical grade foam that protects the stoma from dirt and dust entering the airway, whilst warming and humidifying the incoming air.



Order ref: CR3853 - 65 x 65mm   
Pack size: 30

### Voice Prosthesis Brushes

#### - Standard

Voice Prosthesis Brushes feature soft nylon bristles that prevent damaging the silicone inner surface of the prosthesis. Brushes are tapered to provide effective cleaning of mucus and secretion build-up in the inner lumen of the voice prosthesis.

Order ref: PB/50   
Pack size: 10

## Adhesives and Adhesive Removers

### Acubond Strips

These extra sticky, foam adhesive strips secure the baseplate, giving the ultimate confidence by reducing air leakage.



SB2

Order ref: SB2  
Pack size: 30



### FAHL OPTICLEAR® Adhesive Remover

For laryngectomized patients when using baseplates. Removes adhesives of all types.

Order ref: 33500  
Pack size: 30



### Silbond Silicone Glue - 30ml

Adheres baseplates securely even in the presence of perspiration and moisture.



Order ref: SB1  
Pack size: 1



### Adhere® Silicone Adhesive

Skin friendly, high grade liquid silicone adhesive. Achieves secure attachment of stoma appliances.



Order ref: B-530 - 30ml



### Provox® Luna™ Adhesive

Order ref: 8014  
Pack size: 15



### Provox® Adhesive Strips

Order ref: 8015  
Pack size: 10



### Provox® Glue

Order ref: 7720 - 40ml  
Pack size: 1



### Provox® Adhesive Remover Wipes

Order ref: 8012  
Pack size: 50



### SpiroTect™ Baseplate Adhesive Remover Wipes

These non-sting adhesive removers are designed to give quick and pain-free removal of skin adhesives without causing stinging or skin irritation.



Order ref: PTA600  
Pack size: 30





## Barrier Creams and Lubricants

### Barrier Creams

#### Sorbaderm Barrier Cream



Order ref: 3026 - Sachet - 2g   
Pack size: 20

Order ref: 3028 - Tube - 92g   
Pack size: 1

#### Cutimed PROTECT

*Now available  
in a sachet.*



Order ref: 72652-08 - Sachet - 2g   
Pack size: 20

Order ref: 72652-00 - Tube - 28g   
Pack size: 1

Order ref: 72652-03 - Tube - 90g   
Pack size: 1

### Sterile Lubricating Jelly

#### OptiLube Sterile Lubricating Jelly

The OptiLube sterile lubricating jelly is available in a sachet or tube. The gel is a water based, water soluble gel, which has a non-sticky formulation and is safe and easy to use.



Order ref: 1120 - Sachet - 5g   
Pack size: 1

Order ref: 1121 - Tube - 42g   
Pack size: 1

Order ref: 1122 - Tube - 82g   
Pack size: 1

## Barrier Film Wipes

These wipes provide a protective layer to the skin, and help to keep your skin free from irritation and soreness.

### FAHL OPTIGARD® Protective Dressing Wipes

Preparative skin care for laryngectomized patients when using baseplates or skin adhesives of all types.

Order ref: 33600  
Pack size: 30



### Kapitex® Trachi-Wipe No-Sting Wipe

Trachi-Wipe skin wipes have been specially formulated for neck breathers to help protect the skin around the stoma and tube site from secretions.



Order ref: LA NNZ 3001  
Pack size: 30



### Sorbaderm No-Sting Barrier Film

- Provides up to 72-hours skin protection
- Alcohol free (does not sting)
- Can be used on both intact and damaged skin
- Transparent to allow for continuous visualisation of skin area
- Foam Applicator

Order ref: 3021 - 1ml  
Pack size: 5



Order ref: 3022 - 3ml  
Pack size: 5



### SpiroTect™ Baseplate Barrier Wipe

- Contains a special non-sting formula, does not contain alcohol
- Leaves a protective and quick drying layer on the skin
- Made of superior quality tissue which is softer and stronger
- Helps to reduce skin irritation from baseplates



Order ref: PTA500  
Pack size: 30



### Provox Skin Barrier

Order ref: 8011  
Pack size: ???????



## Voice Prosthesis Brushes

### Blom-Singer® Indwelling Cleaning Brushes

Brushes are designed to fit 16 Fr. and 20 Fr. indwelling voice prostheses ONLY.

Order ref: IN4064 Pack Size: 3	4 - 8mm	<input type="checkbox"/>
Order ref: IN4065 Pack Size: 3	9 - 16mm	<input type="checkbox"/>
Order ref: IN4066 Pack Size: 3	18 - 22mm	<input type="checkbox"/>



### The Laryvox® Duo Brush

Intended for daily cleaning of the voice prosthesis: includes Laryvox® Clip to allow the fixation of the Laryvox® Stoma Light.

Order ref: 29110-06 Pack Size: 6	4 - 6mm	<input type="checkbox"/>
Order ref: 29110-08 Pack Size: 6	7 - 8mm	<input type="checkbox"/>
Order ref: 29110-10 Pack Size: 6	9 - 10mm	<input type="checkbox"/>
Order ref: 29110-12 Pack Size: 6	11 - 13mm	<input type="checkbox"/>
Order ref: 29110-18 Pack Size: 6	14 - 18mm	<input type="checkbox"/>



The Laryvox® Duo Brush comes in packs of 6

### Provox® XL Brushes

Order ref: 7225 Pack size: 6	<input type="checkbox"/>
---------------------------------	--------------------------

### Voiceline Voice Brushes

Order ref: LA PCB 0003 Pack size: 10	<input type="checkbox"/>
---	--------------------------


### Provox® Brushes

Order ref: 7204 Pack size: 6	<input type="checkbox"/>
---------------------------------	--------------------------

## Tube Cleaning Brushes and Swabs

### Kapitex® Tracheostomy Tube Cleaning Brushes

Kapitex® cleaning brushes ensure that the outer and inner tubes are kept free of dry mucus and debris. The brushes are produced with nylon bristles, spirally bound in stainless steel and moulded into a plastic handle. Brush head length is 90mm. NB: These brushes are not suitable for tracheo PVC or silicone tubes.



### Kapitex® Trachi-Swabs

#### - Small - 3.5 x 5mm tubes

Order ref: TR SWB 0150   
Pack size: 50

#### - Large - 5 x 7mm tubes

Order ref: TR SWB 0250   
Pack size: 50

#### - XLarge - 8mm + tubes

Order ref: TR SWB 3050   
Pack size: 50



#### - Orange - 8mm O.D.

Order ref: TR ACC 0001   
Pack size: 2

#### - Blue - 10mm O.D.

Order ref: TR ACC 0002   
Pack size: 2

#### - Purple - 12mm O.D.

Order ref: TR ACC 0003   
Pack size: 2

#### - Green - 14mm O.D.

Order ref: TR ACC 0004   
Pack size: 2

### SpiroTect™ Tube Cleaning Swabs

- Soft foam material prevents damage to tube
- Long and flexible handle allows for a quick and easy clean
- Disposable and single use
- Unlike normal cotton swabs, will not leave loose fibres in the tube

Order ref: PTA300   
Pack size: 50



## Dressings

There are many different types of dressings available that help reduce the risk of infection.

### ActivHeal® Foam Tracheostomy Dressing

The soft, flexible properties of ActivHeal® non-adhesive foam allow the dressing to be worn underneath clothing comfortably and discreetly.

#### ActivHeal® tracheostomy foam square - 10 x 10cm

Order ref: 10009118



### AdvaZorb `T` Lite and Plus `T` Dressings

The AdvaZorb `T` Light is perfect for the treatment of light/medium exuding wounds, whilst the Plus `T` has improved strength and absorbance for high/medium exuding wounds.

#### Advazorb T Lite - 8 x 8cm

Order ref: CR/4212



#### Advazorb T Plus

Order ref: CR/3900



### Lyof foam® Max

Soft and comfortable, Lyof foam Max is designed for moderate to high exuding wounds. Available in a variety of sizes.

#### Lyof foam Max T Square - 9 x 9cm

Order ref: 603207



### Metalline

- Aluminised wound pad
- Flexible and soft
- Good absorbtion
- Permeable to air and water vapour

### Trachi-Dress

Absorbs mucus, reducing the risk of mucus sitting on the skin underneath the tube. The soft, round edges are gentle on the skin and the keyhole is pre-cut for easy fitting around the tracheostomy tube shaft.

### Trachi-Dress F

Is a high performance polyurethane foam dressing which has been specifically formulated with specific characteristics to encourage an optimal healing environment.

## Dressings

### Metalline

- 80 x 90mm

Order ref: 23094   
Pack size: 50

### TrachiDress-F Dressing - Large

- 100 x 82 X 5mm

Order ref: TR DRF 0002   
Pack size: 10

### Trachi-Dress - Small

- 60 x 82mm

Order ref: TR DRE 0001   
Pack size: 20

### Trachi-Dress - Large

- 80 x 100mm

Order ref: TR DRE 0002   
Pack size: 20

### Dressing Pack Spec 10

Order ref: SD9006   
Pack size: 12

### Irripods

- 25 x 20ml

Order ref: 4002   
Pack size: 25

### Normasol

- 25 x 25ml

Order ref: NORMASOL   
Pack size: 25

### FAHL® Tracheal Compress

- Tracheal compresses with aluminium coating to reduce adhesion of the compress to skin around the tracheostoma
- SLIT variant with slit for easy exchange while the tracheostomy tube is in situ



### FAHL® Tracheal Compress

- 80 x 100mm - with slit

Order ref: 30030   
Pack size: 10

### Non woven swabs

Order ref: 3900   
Pack size: 100

### Non woven swabs

Order ref: SWA264T   
Pack size: 100

## Dressings

### The Sensotrach® Duo Ped Slit

- High-quality tracheal compresses with anatomical rounded shape for maximum comfort
- DUO variants with two layers, for patients with moderate to average production of secretion: the layer close to the skin – marked by a green thread - effectively absorbs secretions
- SLIT variants with slit for exchange while the tracheostomy tube is in situ
- PED variants - for paediatric care

30623



### Sensotrach® Duo Ped Slit

- 60 x 70mm

Order ref: 30623



Pack size: 1

### The Sensotrach® Uno Ped Slit

- High-quality tracheal compresses with anatomical rounded shape for maximum wearing comfort
- UNO variants single-layered, for patients with little production of secretion
- SLIT variants with slit for exchange while the tracheostomy tube is in situ
- Special smaller variant for paediatric care (PED)

30523



### Sensotrach® Uno Ped Slit

- 60 x 70mm

Order ref: 30523



Pack size: 1

### Sensotrach® Duo Tracheal Compress

- 90 x 98mm - with slit

Order ref: 30609



Pack size: 1

### Sensotrach® Uno Tracheal Compress

- 90 x 98mm - with slit

Order ref: 30521



Pack size: 1

## Protectors



The Buchanan range of laryngectomy/tracheostomy protectors have been introduced to offer a wider choice. The stoma protectors give the wearer confidence, not only in 'looking good', but also in its airway protection performance.

With a lightweight, soft to the touch construction, you can be assured of day long comfort. Recommended to be washed no more than 3 times.

### Buchanan® Draw-String Protector

- Small

Order ref: LA BUP 0002 - White   
Pack size: 10

### Buchanan® Draw-String Protector

- Large

Order ref: LA BUP 0001 - White   
Pack size: 10

### Buchanan-Lite® Draw-String Protector

- Small

Order ref: LA BLT 2001 - White   
Pack size: 10

### Buchanan-Lite® Draw-String Protector

- Large

Order ref: LA BLT 1001 - White   
Pack size: 10

### Buchanan-Lite® Draw-String Protector

- Small

Order ref: LA BLT 2003 - Beige   
Pack size: 10

### Buchanan-Lite® Draw-String Protector

- Large

Order ref: LA BLT 1003 - Beige   
Pack size: 10

### Buchanan-Lite® Draw-String Protector

- Small

Order ref: LA BLT 2002 - Blue   
Pack size: 10

### Buchanan-Lite® Draw-String Protector

- Large

Order ref: LA BLT 1002 - Blue   
Pack size: 10



## Protectors

Available in a choice of two sizes, large and small. Specifically designed to protect the airway of laryngectomy and tracheostomy patients and to improve their breathing.

### SpiroTect™ Protectors - Antibacterial

- Small - PTA700
- Large - PTA702

- Contain antibacterial treated fabric which can help to reduce infection and filter the inhaled air
- Protect the trachea from inhalation of debris, dust particles and airborne microbes
- Suitable for day and night use
- Made from skin friendly polyamide material with cotton polyester ties for added comfort
- Size: Large 22 x 19cm  
Small 14 x 11cm
- Hand wash up to 5 times



PTA700

Order ref: PTA700	<input type="checkbox"/>
Pack size: 10	
Order ref: PTA702	<input type="checkbox"/>
Pack size: 10	

## Protectors - Velcro Fastening



LA DNX 0001

### Buchanan® DeltaNex Protector

The design of the DeltaNex allows you to attach the protector to a range of attractive cravats. Each protector measures 125mm (w) x 125mm (h). Pack includes 2 neckbands.

### Buchanan® Deltanex Protector

Order ref: LA DNX 0001 - White   
Pack size: 10

### Hirst Bib Protector

The Hirst Bib Protector is made from two layers of a knitted net material with an attached neckband and velcro type fastening.

### Hirst Bib Protector

Order ref: TRA265V - White   
Pack size: 10

### Nex T Collar

- Protects and prevents dust from entering the stoma
- Filter fabric layer aids humidification
- Poly/Acrylic fabric
- Adjustable hook and loop closure fits up to 16.5" (42cm) neck size
- Available in three colours

### Nex T Collar

Order ref: HNLTB - Black   
Pack size: 2

### Nex T Collar

Order ref: HNLTW - White   
Pack size: 2

### Nex T Collar

Order ref: HNLTG - Grey   
Pack size: 2

## Shower and Foam Protectors



Laryngofoam disposable filters warm, moisten and filter the incoming air and can be applied in seconds.

### Laryngofoam

- Beige - 51 x 65mm

Order ref: LA LFF 0003   
Pack size: 30

### Laryngofoam

- Beige - 65 x 65mm

Order ref: LA LFF 0023   
Pack size: 30

### Laryngofoam

- White - 51 x 65mm

Order ref: LA LFF 0001   
Pack size: 30

### SpiroTect™ Universal Shower Protector



PTA101 Ideal for non HME users

Order ref: PTA101   
Pack size: 2

### Provox® Luna™ Shower Aid

Order ref: 7260   
Pack size: 1

### SpiroTect™ Shower Protector

Use with PTC1004 Touch2Speak Cassette, Antibacterial. See page 18.



PTA100

Shower protector with HME in place (HME - PTC1004)

The SpiroTect™ range of easy to use shower protectors enable you to shower without water entering the trachea, taking the hassle out of showering with a stoma.

Order ref: PTA100   
Pack size: 2

## Trachi-Naze +Plus Restoration System

### Filter Cassettes

Each filter cassette contains three filter layers which when sealed onto the skin create a 'closed' system ensuring that all air is warmed, humidified and filtered, restoring a functional airway resistance.

### The pre-filter layer

Filters out dust and larger airborne particles.

### The activated carbon layer

Helps reduce bacteria entering the airway and further cleanses the air of smoke, gases and small particles.

### The heat moisture exchange layer

Provides the heat and moisture needed to keep the tracheal mucosa in good condition. When placed together and in combination, with the system sealed onto the skin, these materials provide the functional airway resistance which is vital for efficient breathing.

### Blue Filter (Night)

Provides normal airway resistance which improves oxygenation. These filters can be worn at night and are easily tolerated during sleep.

### Green Filter (Day)

For normal activity. It provides half the level of airway resistance, but using with the Blue and Green filter, maximises respiratory performance.



Order ref: LATNP 1001   
Pack size: 30



Order ref: LATNP 1002   
Pack size: 30

## Tube and Collar Holders

**Marpac** Tube Holders are designed to hold the tracheostomy tube in position. Each holder incorporates a cough reflex section which allows movement whilst staying in position. The soft foam and cotton material ensures continued comfort.



### Marpac

1 Piece Holder - 4-7 Inch Neck

Order ref: 106   
Pack size: 25

### Marpac

1 Piece Holder - 6-12 Inch Neck

Order ref: 107   
Pack size: 25

### Marpac

2 Piece Collar - Adult

Order ref: 201   
Pack size: 10

### Marpac

2 Piece Collar - Neo-Natal/Paediatric

Order ref: 203   
Pack size: 10

**Portex** Velcro® Trachi Collar Holder is specifically designed to hold tracheostomy tubes firmly in place. Soft, strong and lightweight fabric covered foam for total comfort.

### Portex Velcro® Trachi Collar Holder

- Adult

Order ref: 520000   
Pack size: 1

- Paediatric

Order ref: 520001   
Pack size: 1

- Infant

Order ref: 520002   
Pack size: 1

## Tube and Collar Holders

**SpiroTect™** Tube Holders are specially designed for use by recumbent patients and at night. It features two hook and loop, individually adjustable closures, enabling the Laryngectomee to adjust to a comfortable fit.

### SpiroTect™ Tube Holder

- 1 Piece



PTA200

- Easily adjustable with simple to use hook and loop closures
- Quick and simple closure release in case of emergency
- Narrow foam at the front for added discretion

Order ref: PTA200

Pack size: 10



### SpiroTect™ Tube Holder

- 2 Piece - Adult - PTA201



PTA201

- Two-piece for additional secure positioning
- Easily adjustable with simple to use hook and loop closures
- Faster removal in case of emergency provides extra confidence

Order ref: PTA201

Pack size: 10



### SpiroTect™ Tube Holder

- 2 Piece - Child/Small Adult - PTA301

- Made from soft, skin-friendly material to avoid any discomfort or irritation to the neck area
- Easily adjustable to suit each patient with simple to use hook and loop closures
- Closure can be quickly and simply released in case of emergency
- Two piece for additional secure positioning



Order ref: PTA301

Pack size: 10



## Tube and Collar Holders

### Trachi-Hold Tube Holders

Trachi-Hold tube holders provide a secure method of holding a tracheostomy tube in position. Made from a foam/nylon material to prevent soreness and discomfort. The holders use a velcro style attachment.



#### Trachi-Hold Tube Holder

##### - Mini

Order ref: TR ACC 0030   
Pack size: 10

##### - Small

Order ref: TR ACC 0014   
Pack size: 10

##### - Large

Order ref: TR ACC 0012   
Pack size: 10

### Other

#### Nex Two Piece Tube Holder

Order ref: HNTTHA   
Pack size: 1


#### Provox® LaryClip

Order ref: 7669   
Pack size: 8

## Tube Holders

### Optiflausch® K

- Tracheostomy Tube Holder with Velcro closure (K)
- Assures the secure fixation of tracheostomy tubes
- Foam core surrounded by a soft, non-woven fabric
- Two-part tube holder is pleasantly soft & flexible. It can be adjusted to different neck sizes
- The supporting band is approximately 4 mm thick and 25mm wide




Order ref: 32550  
Pack size: 1



### Optiflausch® K Paediatric

- Tracheostomy tube holder for paediatric use, with adherent fastener (K)
- Assures the secure fixation of tracheostomy tubes
- Broad and soft for comfort
- Approx. 4 mm thick and 25mm wide
- The tube holder is easily and quickly adjusted to the desired neck size (2 parts, adjustable).



Order ref: 32555  
Pack size: 1



### NeckFix® Tube Holder - Duo K

- Tracheostomy tube holder - with adherent fastener - white ("Duo K")
- Assures the secure fixation of tracheostomy tubes
- Foam core surrounded by soft non-woven fabric
- Two-part tube holder is pleasantly soft and flexible. It can be adjusted to different neck sizes
- Individually packed



Order ref: 32630  
Pack size: 1





## Other Accessories

### Stoma Studs

Providing a gentle seal within the stoma, these Stoma Studs are perfect for patients who are unable to tolerate skin adhesives, or who are unable to use baseplates due to recessed or uneven stomas.

The Stoma Studs are available in a variety of sizes ensuring gentle but effective stoma security.



#### Short Length Stud 14.5mm

Code:	ID mm	OD mm	
LATNP 2001	8	10.5	<input type="checkbox"/>
LATNP 2002	10	12.5	<input type="checkbox"/>
LATNP 2003	12	14.5	<input type="checkbox"/>
LATNP 2004	14	16.5	<input type="checkbox"/>
LATNP 2005	16	18.5	<input type="checkbox"/>

**Pack Size: 2**

#### Long Length Stud 22.5mm

Code:	ID mm	OD mm	
LATNP 3001	8	10.5	<input type="checkbox"/>
LATNP 3002	10	12.5	<input type="checkbox"/>
LATNP 3003	12	14.5	<input type="checkbox"/>
LATNP 3004	14	16.5	<input type="checkbox"/>
LATNP 3005	16	18.5	<input type="checkbox"/>

**Pack Size: 2**

### Kapi-Gel® Tracheostoma Spacer

Ideal for difficult stoma anatomy, the Kapi-Gel® Tracheostoma spacer provides cushioning behind the laryngectomy tube.

- Acts as a cushion behind the laryngectomy tube flange/neckplate
- Due to its transparency, Kapi-Gel® is not noticeable around the stoma
- Promotes skin care around the stoma and helps reduce soreness around the stoma and tube site
- Kapi-Gel® is a tacky laminated silicone gel available in two thicknesses; 3mm and 5mm and in two hole dimensions 8mm and 12mm



#### Kapi-Gel Spacers - Thin

Order ref: LATNG 2008 - 3 x 8mm   
Pack size: 10

Order ref: LATNG 2012 - 3 x 12mm   
Pack size: 10

#### Kapi-Gel Spacers - Thick

Order ref: LATNG 4008 - 5 x 8mm   
Pack size: 10

Order ref: LATNG 4012 - 5 x 12mm   
Pack size: 10

## Other Accessories

### ProTrach® DualCare™ HME System

Combining a speaking valve with the strength of an HME with a simple *twist*.

The ProTrach® DualCare™ combines a speaking valve with an HME, enabling speaking and support for your lungs with important heating and moisturising.

The switch between the speaking mode and HME mode is extremely quick and easy, with a simple twist of the lid of the speaking valve.

Speaking mode is made by twisting the lid of the speaking valve until it clicks into speaking mode position.

HME mode is made by twisting the lid of the speaking valve until it clicks into HME position. The HME mode allows you to breathe in and out through the HME.

### ProTrach® DualCare HME

Order ref: 7742



Pack size: 30

### ProTrach® DualCare Speaking Valve (Tinted)

Order ref: 7744



Pack size: 1

### ProTrach® DualCare Speaking Valve (Blue)

Order ref: 7755



Pack size: 1

### Provox® TrachPhone HME



7707

- Speech valve function
- Heat and Moisture Exchanger (HME)
- Integrated oxygen connector
- Lightweight, only 2.9g, minimal stress on tracheostomy or tracheal tube
- Small and compact design, with integrated suction port
- In case of coughing or blockage, suction port acts as a relief valve
- Fits on standard 15mm ISO connector

Order ref: 7707



Pack size: 30

## Other Accessories

### Combiphon® O2 Speaking Valve

- Multi-adaptor fits both 15mm and 22mm connectors
- Speaking valve membrane is “bias-closed” with the advantage that no excess exhalation pressure is required to close the valve, resulting in less fatigue and better speaking quality
- Low-profile design
- Individual sterile packaging
- Single patient product

**Combiphon® is only to be used with Tracheostomised patients, with an intact larynx.**

The use on laryngectomised patients is contraindicated.

Order ref: 27131



Pack size: 1

# For Healthcare Professionals

**Special note:**

If you discover Fittleworth do not stock an item but you believe it to be an Appliance that is available on Prescription, please speak with your Fittleworth Customer Service Advisor who will make enquiries as to whether we can obtain it for your patient.



## FAHL Laryngectomy Tubes

### Laryngotec® HME Compatible Silicone Tubes - Standard (Kombi) and Kombi Lingo (fenestrated)

The Laryngotec® Kombi Lingo version has numerous tiled fenestrations. Designed to reduce the potential for tracheal irritation and granulation tissue formation and facilitates voice when used in conjunction with a voice prosthesis.



Fully compatible with a HME cassette or Hands-Free valve system.

### Laryngotec® HME Compatible Silicone Tubes - Standard (Kombi) and Kombi Lingo (fenestrated)

Size:	O.D. neck Flange mm:	I.D. mm:	Outer curve length mm:	Kombi Order Ref:	Kombi Lingo Order Ref:
7/44	10.2	7.2	44	14920-0744 <input type="checkbox"/>	14720-0744 <input type="checkbox"/>
7/62	10.2	7.2	62	14920-0762 <input type="checkbox"/>	14720-0762 <input type="checkbox"/>
7/65	10.2	7.2	65	14920-0765 <input type="checkbox"/>	14720-0765 <input type="checkbox"/>
8/44	11.2	8.2	44	14920-0844 <input type="checkbox"/>	14720-0844 <input type="checkbox"/>
8/62	11.2	8.2	62	14920-0862 <input type="checkbox"/>	14720-0862 <input type="checkbox"/>
8/70	11.2	8.2	70	14920-0870 <input type="checkbox"/>	14720-0870 <input type="checkbox"/>
9/44	12.2	9.4	44	14920-0944 <input type="checkbox"/>	14720-0944 <input type="checkbox"/>
9/62	12.4	9.4	62	14920-0962 <input type="checkbox"/>	14720-0962 <input type="checkbox"/>
9/75	12.4	9.4	75	14920-0975 <input type="checkbox"/>	14720-0975 <input type="checkbox"/>
10/44	13.4	10.4	44	14920-1044 <input type="checkbox"/>	14720-1044 <input type="checkbox"/>
10/62	13.4	10.4	62	14920-1062 <input type="checkbox"/>	14720-1062 <input type="checkbox"/>
10/85	13.4	10.4	85	14920-1085 <input type="checkbox"/>	14720-1085 <input type="checkbox"/>
11/44	14.4	11.4	44	14920-1144 <input type="checkbox"/>	14720-1144 <input type="checkbox"/>
11/62	14.4	11.4	62	14920-1162 <input type="checkbox"/>	14720-1162 <input type="checkbox"/>
12/44	15.4	12.4	44	14920-1244 <input type="checkbox"/>	14720-1244 <input type="checkbox"/>

**Pack Size: 1**

## Kapitex Trachi-Naze® Plus Laryngectomy Tubes

- Designed to hold the Trachi-Naze Plus filter cassettes
- Produced from medical grade silicone
- Soft, supportive and easy to use
- Restores essential nasal airway function

Code:	Cannula Length:	Inner Diameter:	Outer Diameter:	Pack Size:	
LA TNL 0840	40mm	8.0mm	10.5mm	1	☐
LA TNL 1020	20mm	10.0mm	12.5mm	1	☐
LA TNL 1030	30mm	10.0mm	12.5mm	1	☐
LA TNL 1040	40mm	10.0mm	12.5mm	1	☐
LA TNL 1220	20mm	12.0mm	14.5mm	1	☐
LA TNL 1230	30mm	12.0mm	14.5mm	1	☐
LA TNL 1240	40mm	12.0mm	14.5mm	1	☐
LA TNL 1250	50mm	12.0mm	14.5mm	1	☐
LA TNL 1260	60mm	12.0mm	14.5mm	1	☐
LA TNL 1420	20mm	14.0mm	16.5mm	1	☐
LA TNL 1440	40mm	14.0mm	16.5mm	1	☐

## SpiroTect™ Laryngectomy Tubes

### SpiroTect™ Laryngectomy Tube Perforated - PTT2001-2008

- Perforated tube is designed to allow for speech
- Medical grade silicone material
- PVC and Latex free
- SpiroTect™ 44mm tube is the same length as alternative brands 36mm and 62mm is the same length as alternative brands 55mm

Code:	Size:	Length:	Inner Diameter:	Outer Diameter:	Pack Size:
PTT2001	8	44mm	9.5mm	12.0mm	1 <input type="checkbox"/>
PTT2002	9	44mm	10.5mm	13.5mm	1 <input type="checkbox"/>
PTT2003	10	44mm	12mm	15.0mm	1 <input type="checkbox"/>
PTT2004	12	44mm	13.5mm	17.0mm	1 <input type="checkbox"/>
PTT2005	8	62mm	9.5mm	12.0mm	1 <input type="checkbox"/>
PTT2006	9	62mm	10.5mm	13.5mm	1 <input type="checkbox"/>
PTT2007	10	62mm	12.0mm	15.0mm	1 <input type="checkbox"/>
PTT2008	12	62mm	13.5mm	17.0mm	1 <input type="checkbox"/>

### SpiroTect™ Laryngectomy Tube Non-Perforated - PTT2009-2016

Code:	Size:	Length:	Inner Diameter:	Outer Diameter:	Pack Size:
PTT2009	8	44mm	9.5mm	12.0mm	1 <input type="checkbox"/>
PTT2010	9	44mm	10.5mm	13.5mm	1 <input type="checkbox"/>
PTT2011	10	44mm	12mm	15.0mm	1 <input type="checkbox"/>
PTT2012	12	44mm	13.5mm	17.0mm	1 <input type="checkbox"/>
PTT2013	8	62mm	9.5mm	12.0mm	1 <input type="checkbox"/>
PTT2014	9	62mm	10.5mm	13.5mm	1 <input type="checkbox"/>
PTT2015	10	62mm	12.0mm	15.0mm	1 <input type="checkbox"/>
PTT2016	12	62mm	13.5mm	17.0mm	1 <input type="checkbox"/>

## Shiley™ Tracheostomy Tubes with Disposable Inner Cannula

### Disposable Inner Cannula Tubes for Safe and Convenient Tracheostomy Care

Offering safe and convenient tracheostomy care, these disposable inner cannula tubes provide a quick and convenient restoration of a patient's airway.

Description:	Code:	Size:	I.D:	O.D:	Length:	
Tracheostomy tube cuffed with disposable inner cannula - DCT	4DCT	4	5.0mm	9.4mm	62mm	<input type="checkbox"/>
	6DCT	6	6.4mm	10.8mm	74mm	<input type="checkbox"/>
	8DCT	8	7.6mm	12.2mm	79mm	<input type="checkbox"/>
	10DCT	10	8.9mm	13.8mm	79mm	<input type="checkbox"/>
Tracheostomy tube cuffed with disposable inner cannula - Fenestrated DFEN	4DFEN	4	5.0mm	9.4mm	62mm	<input type="checkbox"/>
	6DFEN	6	6.4mm	10.8mm	74mm	<input type="checkbox"/>
	8DFEN	8	7.6mm	12.2mm	79mm	<input type="checkbox"/>
	10DFEN	10	8.9mm	13.8mm	79mm	<input type="checkbox"/>
Tracheostomy tube cuffless with disposable inner cannula - DCFS	4DCFS	4	5.0mm	9.4mm	62mm	<input type="checkbox"/>
	6DCFS	6	6.4mm	10.8mm	74mm	<input type="checkbox"/>
	8DCFS	8	7.6mm	12.2mm	79mm	<input type="checkbox"/>
	10DCFS	10	8.9mm	13.8mm	79mm	<input type="checkbox"/>
Tracheostomy tube cuffless with disposable inner cannula - Fenestrated - DCFN	4DCFN	4	5.0mm	9.4mm	62mm	<input type="checkbox"/>
	6DCFN	6	6.4mm	10.8mm	74mm	<input type="checkbox"/>
	8DCFN	8	7.6mm	12.2mm	79mm	<input type="checkbox"/>
	10DCFN	10	8.9mm	13.8mm	79mm	<input type="checkbox"/>

**Pack Size: 1**



## Shiley™ Reusable Inner Cannula Tracheostomy Tubes

### An economical alternative for long-term and home care patients.

The unique inner cannula can be cleaned an unlimited number of times within a 29-day tube-use cycle. The high-volume, low-pressure cuff provides a gentle and effective seal, reducing the risk of tracheal damage.

Description:	Code:	Size:	I.D:	O.D:	Length:	
Tracheostomy tube cuffed with inner cannula - LPC	4LPC	4	5.0mm	9.4mm	65mm	<input type="checkbox"/>
	6LPC	6	6.4mm	10.8mm	76mm	<input type="checkbox"/>
	10LPC	10	8.9mm	13.8mm	81mm	<input type="checkbox"/>
Tracheostomy tube cuffed with inner cannula - Fenestrated - FEN	4FEN	4	5.0mm	9.4mm	65mm	<input type="checkbox"/>
	6FEN	6	6.4mm	10.8mm	76mm	<input type="checkbox"/>
	8FEN	8	7.6mm	12.2mm	81mm	<input type="checkbox"/>
	10FEN	10	8.9mm	13.8mm	81mm	<input type="checkbox"/>
Tracheostomy tube cuffless with inner cannula - Fenestrated - CFN	4CFN	4	5.0mm	9.4mm	65mm	<input type="checkbox"/>
	6CFN	6	6.4mm	10.8mm	76mm	<input type="checkbox"/>
	10CFN	10	8.9mm	13.8mm	81mm	<input type="checkbox"/>
Tracheostomy tube cuffless with disposable inner cannula - CFS	4CFS	4	5.0mm	9.4mm	65mm	<input type="checkbox"/>
	6CFS	6	6.4mm	10.8mm	76mm	<input type="checkbox"/>
	8CFS	8	7.6mm	12.2mm	81mm	<input type="checkbox"/>
	10CFS	10	8.9mm	13.8mm	81mm	<input type="checkbox"/>

### Pack Size: 1

## Shiley™ Paediatric Tracheostomy Tubes Cuffed and Cuffless

The Paediatric Tracheostomy Tubes provide patients with the right fit, whilst maintaining airway patency, patient safety and comfort.

Description:	Code:	Size:	I.D:	O.D:	Length:	
Paediatric Tracheostomy tube, cuffless - PED	3.0 PED	3.0	3.0mm	4.5mm	39mm	<input type="checkbox"/>
	3.5 PED	3.5	3.5mm	5.2mm	40mm	<input type="checkbox"/>
	4.0 PED	4.0	4.0mm	5.9mm	41mm	<input type="checkbox"/>
	4.5 PED	4.5	4.5mm	6.5mm	42mm	<input type="checkbox"/>
	5.0 PED	5.0	5.0mm	7.1mm	44mm	<input type="checkbox"/>
	5.5 PED	5.5	5.5mm	7.7mm	46mm	<input type="checkbox"/>
Tracheostomy tube, long - cuffless - PDL	5.0 PDL	5.0	5.0mm	7.1mm	50mm	<input type="checkbox"/>
	5.5 PDL	5.5	5.5mm	7.7mm	52mm	<input type="checkbox"/>
	6.0 PDL	6.0	6.0mm	8.3mm	54mm	<input type="checkbox"/>
	6.5 PDL	6.5	6.5mm	9.0mm	56mm	<input type="checkbox"/>
Tracheostomy tube, cuffed - PDC	4.0 PDC	4.0	4.0mm	5.9mm	41mm	<input type="checkbox"/>
	4.5 PDC	4.5	4.5mm	6.5mm	42mm	<input type="checkbox"/>
	5.0 PDC	5.0	5.0mm	7.1mm	44mm	<input type="checkbox"/>
	5.5 PDC	5.5	5.5mm	7.7mm	46mm	<input type="checkbox"/>
Tracheostomy tube, long - cuffed - PLC	5.0 PLC	5.0	5.0mm	7.1mm	50mm	<input type="checkbox"/>
	5.5 PLC	5.5	5.5mm	7.7mm	52mm	<input type="checkbox"/>
	6.0 PLC	6.0	6.0mm	8.3mm	54mm	<input type="checkbox"/>

**Pack Size: 1**

## Jaw Rehabilitation

TheraBite® is a portable device that provides patients with an anatomically correct passive jaw motion. It also helps reduce patient anxiety by allowing them to control the extent and length of each stretch. The TheraBite® system provides both anatomically correct stretching and passive motion therapy for effective jaw rehabilitation. Available in two sizes: Regular for adults and paediatric for smaller mouths.

### Contents:

- 1 x Jaw Mobiliser
- 30 x Range of Motion Scale
- 1 x Patient Progress
- 4 x Bite Pads
- 1 x Carry Bag
- 1 x Hand Aid
- 1 x Instructions for use

Description:	Code:	Size:	Pack Size:
TheraBite® Jaw Motion Rehab System™	TH001	Adult	1 <input type="checkbox"/>
TheraBite® Jaw Motion Rehab System™	TH002	Paediatric	1 <input type="checkbox"/>
TheraBite® Bite Pads	PA001	Regular/Adult	4 <input type="checkbox"/>
TheraBite® Bite Pads	PA002	Endentulous	4 <input type="checkbox"/>
TheraBite® Bite Pads	PA003	Paediatric	4 <input type="checkbox"/>

### Provox Larytube Fenestrated (including 5 HME Cassette HiFlow and 1 Provox Brush)

Code:	Cannula Length:	Inner Diameter:	Outer Diameter:	Pack Size:
7637	36mm	9.5mm	12.0mm	1 <input type="checkbox"/>
7638	55mm	9.5mm	12.0mm	1 <input type="checkbox"/>
7640	36mm	10.5mm	13.5mm	1 <input type="checkbox"/>
7641	55mm	10.5mm	13.5mm	1 <input type="checkbox"/>
7643	36mm	12.0mm	15.0mm	1 <input type="checkbox"/>
7644	55mm	12.0mm	15.0mm	1 <input type="checkbox"/>
7646	36mm	13.5mm	17.0mm	1 <input type="checkbox"/>
7647	55mm	13.5mm	17.0mm	1 <input type="checkbox"/>



# Home Delivery Service

## About Home Delivery - How does this service work?

You can order through Fittleworth in six ways, all of which are very simple and all of which are FREE:

### **NHS Electronic Prescription Service:**

We accept electronic prescriptions (EPS). Your GP can send your prescription electronically to us at Fittleworth. To find out more, please call, text or email us using the details below or contact your GP.

### **By Post:**

Send your prescription in one of the envelopes that we provide or to our address: FREEPOST, Fittleworth.

### **On our Freephone Telephone Numbers:**

National: **0800 378 846**

Scotland: **0800 783 7148**

Nurse Line: **0800 378 413**

### **Text Messaging:**

Text: **07800 005 658**

### **Email:**

You can email your order to us at:

**orders@fittleworth.com**

### **Online:**

**www.fittleworth.com**

## Reordering

To make it easier to reorder each month, the space below can be used to list the items that your patient requires:

<b>Product description</b>	<b>Order ref</b>	<b>Quantity</b>

## Reordering

Product description	Order ref	Quantity



**FREEPOST**  
**Fittleworth**



National:  
**0800 378 846**  
Scotland:  
**0800 783 7148**  
Nurse Line:  
**0800 378 413**  
Text:  
**07800 005 658**



[www.fittleworth.com](http://www.fittleworth.com)



NOVEMBER 2020